

# Three group classification of participants based on fully automated plasma $\beta$ -amyloid measurements to achieve high positive and negative predictive values

**Kazuto Yamashita<sup>1</sup>**, Masahiro Miura<sup>1</sup>, Kota Nagai<sup>2</sup>, David Verbel<sup>3</sup>, Shigeki Iwanaga<sup>1</sup>, Toshiyuki Sato<sup>1</sup>, Tomokazu Yoshida<sup>4</sup> and Atsushi Iwata<sup>5</sup>  
 (1) Central research laboratories, Sysmex Corporation, (2) Japan and Asia Clinical Development Department, Eisai Co., Ltd., (3) Biostatistics, Eisai Inc., (4) Sysmex Corporation (5) Department of Neurology, Tokyo Metropolitan Geriatric Hospital and Institute of gerontology

## Background

- The accelerated development of disease-modifying therapies targeting  $\beta$ -amyloid ( $A\beta$ ) for the treatment of Alzheimer's disease (AD) has increased the need for blood-based biomarkers to predict brain  $A\beta$  status.
- We reported that plasma  $A\beta_{1-42}$  ( $A\beta_{42}$ ) to  $A\beta_{1-40}$  ( $A\beta_{40}$ ) ratio measured by our fully automated immunoassay platform (HISCL™ series) predicted brain  $A\beta$  status defined by amyloid positron emission tomography (PET) as assessed by Centiloids (CL) with area under the curve (AUC) of 0.922.
- In the previous analysis, we determined a cut-off value of 0.102 using the Youden index. Using this cut-off value, we achieved high negative predictive value (NPV) of 94.0% and moderate positive predictive value (PPV) of 79.6%.
- In this study, we classified participants into three groups (plasma  $A\beta$  positive, intermediate, and negative groups) by their plasma  $A\beta_{42}/A\beta_{40}$  ratio, in order to improve PPV of our plasma  $A\beta$  assay.



- Rapid measurement**  
17 minutes per test
- Small sample volume**  
30  $\mu$ L per test
- High throughput**  
100 or 200 tests per hour (depending on the instrument)

Figure 1. Appearance and features of HISCL series.

## Methods

- Plasma  $A\beta_{40}$  and  $A\beta_{42}$  were measured using a fully automated platform in a set of plasma samples sourced from participants in the screening period of the elenbecostat Phase 3 program.
- In this analysis, we included 172 amyloid PET positive and 199 negative participants.
- Brain  $A\beta$  status was determined by predefined Centiloid unit (CL) of amyloid PET (cut-off value defined previously as 32.21 CL).
- The AUC for plasma  $A\beta_{42}/A\beta_{40}$  ratio for predicting brain  $A\beta$  status was obtained by performing ROC analysis using logistic regression.
- We determined the cut-off value of our plasma  $A\beta_{42}/A\beta_{40}$  ratios that would result in a PPV of 90.0% or more. We then utilized this cut-off value and the prior reported cut-off value of 0.102 to divide participants into plasma  $A\beta$  positive, intermediate, and negative groups.




Figure 2. Assay principle of HISCL series.

## Results

### 1. Participants demographics

- Amyloid PET positive group was significantly older, more often carried one or two *APOE*  $\epsilon$ 4 allele, and more often used florbetaben as an amyloid PET probe.

Table 1. Participants demographics

	Amyloid PET negative (n=199)	Amyloid PET positive (n=172)	P value
Age (y), mean $\pm$ SD	69.1 $\pm$ 8.8	72.8 $\pm$ 7.1	< 0.001
Sex, female/male	109/90	91/81	NS
Race, White/Japanese/Other	130/48/21	111/51/10	NS
<i>APOE</i> $\epsilon$ 4 status, -/+ /NA	160/38/1	68/103/1	< 0.001
Amyloid PET probe, FBB/FBP/FMM	138/16/45	100/28/44	0.024
MMSE, mean $\pm$ SD	26.4 $\pm$ 1.7	26.5 $\pm$ 1.9	NS
CDR-SB, mean $\pm$ SD	2.4 $\pm$ 1.0	2.5 $\pm$ 1.0	NS
Clinical disease staging, MCI due to AD/mild AD	40/159	35/172	NS

The differences between groups were evaluated using the Student's t-test for continuous measures and the  $\chi^2$  test for categorical variables.

Abbreviations: y, years; SD, standard deviation; NS, not significant; *APOE*, Apolipoprotein E; NA, not available; FBB, florbetaben; FBP, florbetapir; FMM, flutemetamol; MMSE, Mini-Mental State Examination; CDR-SB, Clinical Dementia Rating-Sum of Boxes; MCI, Mild Cognitive Impairment.

### 2. ROC analysis

- Plasma  $A\beta_{42}/A\beta_{40}$  ratio measured using HISCL series predicted brain  $A\beta$  status defined by amyloid PET with an AUC of 0.93.




Figure 3. Clinical performance of plasma  $A\beta_{42}/A\beta_{40}$  ratio.

The AUC is described with 95% confidence interval (CI) which is calculated based on DeLong method.

### 3. Three group classification

- A cut-off value of 0.092 was determined as a value to achieve a PPV of at least 90.0%.
- PPV in positive  $A\beta$  group was 90.1% while NPV in the negative  $A\beta$  group was 95.8%.




Figure 4. Three group classification using plasma  $A\beta_{42}/A\beta_{40}$  ratio. Participants were classified into positive ( $\leq 0.092$ ), intermediate ( $> 0.092$  and  $\leq 0.102$ ), and negative ( $> 0.102$ )  $A\beta$  groups based on plasma  $A\beta_{42}/A\beta_{40}$  ratio.

Table 2. Predictive values of plasma  $A\beta_{42}/A\beta_{40}$  ratio.

Plasma $A\beta$ group	Amyloid PET		Predictive value	Frequency
	Negative	Positive		
Negative ( $A\beta_{42}/A\beta_{40} > 0.102$ )	158	7	<b>95.8% (NPV)</b>	44.5%
Intermediate ( $0.092 < A\beta_{42}/A\beta_{40} \leq 0.102$ )	30	65	68.4% (PPV)	25.6%
Positive ( $A\beta_{42}/A\beta_{40} \leq 0.092$ )	11	100	<b>90.1% (PPV)</b>	29.9%

## Conclusion

- Three group classification allowed our plasma  $A\beta$  assay to achieve PPV and NPV  $\geq 90\%$  with 74.4% of participants classifiable as  $A\beta$  positive or negative groups.
- This result indicated that our assay may contribute to reduce amyloid PET scan or CSF  $A\beta$  testing, which could be helpful in applications such as the recruitment step of clinical trials.
- However, it should be noted that predictive values and frequency will vary depending on the prevalence of amyloid PET positive participants.